

BioClave 12L

- ✦ N-type semi automatic Autoclave
- ✦ Front loading bioclave with a capacity of 12 Ltr
- ✦ Microprocessor controlled
- ✦ Ensures biosafety & Sterilization of non porous instruments
- ✦ Completes dry cycle with semi - opened door
- ✦ Door lock for safety
- ✦ Stainless steel chamber and 2 stainless steel trays for surgical instruments
- ✦ Silicone ring seals the chamber
- ✦ Thermostat prevents over heating
- ✦ Diameter of the chamber 21.4 cm
- ✦ Dimension - 40.7 cm width x 38.5 high x 46.5 cm deep
- ✦ Net weight 23 Kg

GNATUS

AUTOCLAVE

