

BioClave 21L

- ✦ N-type semi automatic Autoclave
- ✦ Front loading bioclave with a capacity of 21 Ltr
- ✦ Microprocessor controlled
- ✦ Ensures biosafety & Sterilization of non porous instruments
- ✦ Completes dry cycle with semi-opened door
- ✦ Door lock with safety system against accidental doors opening
- ✦ Stainless steel chamber and 3 stainless steel trays for surgical instruments / Implant Kit
- ✦ Silicone ring seals the chamber
- ✦ Thermostat prevents over heating
- ✦ Complete cycle 60 minutes
- ✦ 43.5 cm width x 41.0 high x 58.0 cm deep
- ✦ Net weight 30 Kg

GNATUS

AUTOCLAVE

