



Classe A3 Plus

CLASSE A3 PLUS

- ✦ Option available for extra wide or narrow backrest. Seat with seamless upholstery for total patient comfort with double articulated headrest
- ✦ Movable right arm rest
- ✦ Aluminum alloy chair to prevent rusting
- ✦ Synchronized movement of seat and backrest
- ✦ 7 Programs in the chair which include - 4 programmable working positions, preset auto return to zero, spitting and last working position
- ✦ Trendelenburg position for dental emergencies
- ✦ Independent up/down chair movements along with back/forth backrest movements
- ✦ Safety system on chair base and backrest and assistant side to prevent entrapment
- ✦ Chair base can also be used for operating the suctions
- ✦ Instrument control panel incorporates with 3 figure numerical display for greater data reading precision



CE ISO 9001:2008



Anthos Chair comes with Suction System VS 250 S

VENUS PLUS L LED

- ✦ 3 axis rotation
- ✦ Sensor operated
- ✦ Inhouse temp. : 5,000k
- ✦ Light intensity : 3,000 to 50,000 Lux
- ✦ Power absorption : 10 W

- ✦ Anthos dental stool provided with chair.
- ✦ Spitton bowl can be dismantled completely with removeable cup & bowl fill unit

CONTINENTAL DELIVERY UNIT

- ✦ 5 modules with lockable hangers
- ✦ Separate feather-touch keypad
- ✦ Each instrument has a retractions system
- ✦ X-Ray viewer with automatic turn off
- ✦ Removable & autoclavable resting silicon pads
- ✦ Movable instrument tray below the delivery unit

WATER UNIT AND ASSISTANT ARM

- ✦ High quality ceramic spittoon with programmable water input for spittoon and cup
- ✦ Assistant Instrument table, fitted on articulated arm for 2 suction and optional instruments

