



# Classe A5

Anthos Classe A5

- ✦ Extra wide back rest & seat with seamless upholstery for total patient comfort with double articulates headrest
- ✦ Aluminum alloy chair to prevent rusting
- ✦ Movable right arm rest
- ✦ Synchronized movement of seat & backrest
- ✦ 7 Programs in the chair which include...4 programmable working positions, preset auto return to zero, spitting & last working position
- ✦ Trendelenburg position for dental emergencies
- ✦ Independent up / down chair movements along with back/forth backrest movements
- ✦ Safety system on chair base, backrest & assistant side to prevent entrapment
- ✦ Chair base can also be used for operating the suctions (high & low)
- ✦ The LCD display illustrates instruments speed & other several icons provide information on dental unit / operations.



Venus Plus L LED Light



CE ISO 9001:2008

DÜRR DENTAL

Anthos Chair comes with Suction System VS 250 S



### CONTINENTAL DELIVERY UNIT

- ✦ 5 modules with lockable hangers
- ✦ Separate feather-touch keypad
- ✦ The O.D.R. (Out Drop no Retraction) device is a standard feature which avoids retraction of liquid droplets on instrument tips
- ✦ X-Ray viewer with automatic turn off
- ✦ Removable & autoclavable resting silicon pads
- ✦ Movable instrument tray below the delivery unit
- ✦ Anthos micromotor & scaler
- ✦ 2 nos. mid-west connections for high speed handpieces

### WATER UNIT AND ASSISTANT ARM

- ✦ High quality ceramic spittoon with programmable water input for spittoon & cup
- ✦ Assistant instrument table, fitted on articulated arm for 2 suction (Low & High supplied with DÜRR VS 250 S suction system) & optional instruments

### VENUS PLUS L LED LIGHT

- ✦ The dual-reflector LED operating light belongs to a new generation of intelligent lighting; it consumes up to 10 times less than a traditional light & has a guaranteed working life of over 50,000 hours
- ✦ 3-axis rotation for greater comfort & precision
- ✦ Adjustment of light intensity from 3,000 to 50,000 Lux
- ✦ Easy to clean removable handle

