

# DCI Edge 360



**2** Years  
Warranty

**Provision for right & left  
hand dentistry**



- + Aluminum casting chair
- + New concept & aesthetic design
- + Seamless & soft upholstery for patient comfort
- + Synchronised movement of seat & backrest
- + Chair can lift upto 250 kg weight
- + Hydraulic delivery system from Japan
- + Hand rest - with two direction rotation
- + Safety with protection against any obstruction when chair goes down
- + Double articulated headrest with pull & push button for pediatric & handicapped patients
- + All connection in the base of chair

## **WATER UNIT**

- + 90° Movements best for pedo & old age patients
- + Removable cuspidor, easy to be disinfected
- + Air & water suction filter to prevent drainage block

## **Assistant Side**


- + Three points - one for 3 way syringe & 2 for low high vacuum suction
- + Chair control on assistant side.

## **DENTIST STOOL**


- + With adjustable up & down position & also backrest can be adjusted upto 250 degree

## **DELIVERY UNIT**

- + Delivery unit with 3 air rotor point & one point for 3-way syringe
- + All the control of chair are on unit with touch panel



# DCI Edge



The seat can be rotated about 30 degree, More convenient for maintenance and multiple operations.

- ✦ 3 axis rotation (For easy positioning and projection)
- ✦ 3-intensity : 16000 Lux, 26000 Lux & 38000 Lux
- ✦ Precise illumination Pattern : 3" X 6" (Prevents Patient Glare)
- ✦ Color Temperature : 5000K (Accurate shade matching)
- ✦ Color Rendering Index (CRI) : 90
- ✦ Shadow free light (Less eye fatigue & Increase Visibility)
- ✦ Reflective technology : Uses 500 reflective facet (Uniform light pattern)
- ✦ Durable Design : Heavy-duty structure & handle design
- ✦ Hermetically sealed front screen makes external cleaning easier and stops any dust or spring fluid getting inside the light

# DCI LIGHT

