

- + Simple, easy & quick to use
- + Convenient cordless device using battery
- + Cordless & electric motorized filling system
- + Disposable fill needle
- + 360° rotation at both cartridge needle & pack tip
- + High capacity lithium-ion battery

Clinical Benefit

- + Less time consuming
- + Less microbial leakage
- + Adapts better to grooves & depressions at chronol. Walls & Lateral Canals
- + No dead space so no secondary injection

Pack :

- + 360° rotation without taking out the pack tip
- + Outstanding heating performance (Max Temp. : 300°C)
- + Ergonomically designed body for comfortable grip
- + Highly visible LED display
- + "Hold" switch for safety
- + Automatic off after 5 minutes
- + Dimension : 153mm x 27mm x 25mm (WxLxH)
- + Weight : 90 g (Including battery)
- + Operating Voltage : 3.7V
- + Range of setting temperature : 140~300° C
- + Auto off : 5min

Fill :

- + Easy-Install fill needle
- + High visible LED display
- + 360° rotating cartridge needle
- + Fully motorized obturation for smooth & easy filling
- + Ergonomically designed body for comfortable grip
- + Disposable fill needle
- + Auto reverse function when the cartridge is empty
- + "Hold" switch for safety
- + Dimension : 223mm x 27mm x 26mm (WxLxH)
- + Weight : 177 g (Including battery)
- + Operating Voltage : 3.7V
- + Range of setting temperature : 140~200° C
- + Auto off : 5min

Genuine System, Genius Result

GENESYS

Genuine Root Canal Obturation system for Thermal 3D Obturation of Root Canal



Pack

Fill

POWERED BY

**UNICORN
DenMart** 

Behind The Best Indian Dentists