

A DREAM MACHINE

For Complete Root Canal Treatment

Four-in-One (Endo Motor, Apex Locator, Filler & Obturator System). Complete Solution for an Endodontist

Features

- ✦ One touch color wide screen
- ✦ Compact & light weight
- ✦ Endo motor & Apex Locator at the same time
- ✦ Options to Mount on dental chair or movable cart type
- ✦ Best accuracy & systematic diagnosis

Endo Motor Features

- ✦ Auto-reverse mode
- ✦ Torque control from 0.6 to 6.4Ncm
- ✦ Long shelf life with swiss motor
- ✦ 10 Memory saving function
- ✦ Reciprocating mode

Obturator (Fill) Features

- ✦ Easy install for fill needle
- ✦ Fixed cartridge needle
- ✦ Fully motorized obturation for smooth filling
- ✦ Long shelf-life with swiss motor
- ✦ Comfortable grip with ergonomic design
- ✦ Auto reverse function when cartridge is empty
- ✦ Automatic off after 5 minutes
- ✦ Option to reused cartridge

Apex Locator Features

- ✦ Best accuracy in any root canal
- ✦ Self test for system failure
- ✦ Simple, easy to use
- ✦ Color display

Obturator (Pack) Features

- ✦ 360° rotating for pack tip
- ✦ Outstanding heating performance (Max. 4000C)
- ✦ Comfortable grip
- ✦ Hold switch for safety
- ✦ Temperature range (140~400°C) displayed in LED window
- ✦ Automatic off after 5 minutes



Dental Unit For Endodontic Treatment

4 IN 1 SYSTEM

ENDO MOTOR

Shape the Root Canal

APEX LOCATOR

Measure The Working Length

OBTURATOR PACK

Pack Root Canal

OBTURATOR FILL

Fill Root Canal

POWERED BY

UNICORN DenMart 

Behind The Best Indian Dentists