

2 Years Warranty



Victor V286 Chair comes with VMax 450 Suction System



- + Aluminum alloy chair to prevent rusting & reliability
- + Synchronized movement of seat and backrest
- + 6 Programs in the chair which include - 3 programmable working positions, preset auto return to zero, spitting and last working position
- + Trendelenburg position for dental emergencies
- + Independent up/down chair movements along with back/forth backrest movements
- + Instrument control panel incorporates with 3 figure numerical display for greater data reading precision

LED LIGHT

- + 3 axis rotation
- + Sensor operated
- + Inhouse temp. : 4,500k
- + Light intensity : 30,000 Lux

CONTINENTAL DELIVERY UNIT

- + 5 modules with lockable hangers
- + Separate feather-touch keypad

- + Each instrument has a retractions system
- + X-Ray viewer with automatic turn off
- + Removable & autoclavable resting silicon pads
- + Movable instrument tray below the delivery unit

WATER UNIT AND ASSISTANT ARM

- + High quality ceramic spittoon with programmable water input for spittoon and cup
- + Assistant Instrument table, fitted on articulated arm for 2 suction and optional instruments
- + Dental stool with chair
- + Spittoon bowl can be dismantled completely with removable cup & bowl fill unit

MULTI FUNCTIONAL FOOT CONTROL

- + Joy stick for patient chair control
- + Rinse & Zero position
- + Light & cup fill control
- + Chip air & chip water control
- + Instrument speed & power control