

- ✦ New compact design, new “free stand” unit body
- ✦ New light weight medical arm, compact shape, maneuverable, great mobility with great flexibility to ergonomics
- ✦ Option for continental or side delivery unit
- ✦ Extra wide back rest & seat with seamless upholstery for total patient comfort with double articulates headrest
- ✦ Aluminum alloy chair to prevent rusting
- ✦ 7 Programs in the chair which include...4 programmable working positions, preset auto return to zero, spitting & last working position
- ✦ Trendelemburg position for dental emergencies
- ✦ Safety system on chair base, backrest & assistant side to prevent entrapment
- ✦ The LCD display illustrates instruments speed & other several icons provide information on dental unit / operations.



2 Years Warranty



CE ISO 9001:2008 EN ISO 13485:2012



VENUS PLUS L LED LIGHT

- ✦ The dual-reflector LED operating light belongs to a new generation of intelligent lighting; it consumes up to 10 times less than a traditional light & has a guaranteed working life of over 50,000 hours
- ✦ 3-axis rotation for greater comfort & precision
- ✦ Adjustment of light intensity from 3,000 to 50,000 Lux
- ✦ Easy to clean removable handle

CONTINENTAL DELIVERY UNIT

- ✦ 5 modules with lockable hangers
- ✦ Separate feather-touch keypad
- ✦ The O.D.R. (Out Drop no Retraction) device is a standard feature which avoids retraction of liquid droplets on instrument tips
- ✦ X-Ray viewer with automatic turn off
- ✦ Removable & autoclavable resting silicon pads
- ✦ Movable instrument tray below the delivery unit
- ✦ Anthos micromotor & scaler
- ✦ 2 nos. mid-west connections for high speed handpieces

WATER UNIT AND ASSISTANT ARM

- ✦ High quality ceramic spittoon with programmable water input for spittoon & cup
- ✦ Assistant instrument table, fitted on articulated arm for 2 suction (Low & High supplied with DÜRR VS 250 S suction system) & optional instruments