

RAY68(M)

Floor Mount




X-RAY MACHINE

- ✦ Option to choose between Floor & Wall Mount (A/C model)
- ✦ Soft Scissor arm
- ✦ 70 kVp/7mA X-Ray tube
- ✦ Focal Spot 0.8 mm
- ✦ Focal Length 20 cm
- ✦ Inherent Filtration ≥ 0.5 mm AL
- ✦ Exposure Time 0.06 to 2.0 second
- ✦ Tube Head & cone are internally LEAD coated
- ✦ Remote control for exposure
- ✦ Microprocessor controlled



RAY68(S)

Wall Mount



CE ISO 13485:2012

