

SEA -16N

- + Temperature digital display
- + Imported Italy Sirai Valve
- + Gasket is easy to be cleaned
- + Temperature deviation is less than 0.5° C
- + 1800W heating elements, speed up the heating, more efficient
- + Voltage : 220 V or 120 V
- + Frequency : 50Hz/60Hz
- + Power : 1800W
- + Chamber Capacity : 16L
- + Sterilization temperature : 121° C/20 min, 134° C/6 min
- + Chamber Size : φ200 x 360mm



SEA -45L / 29L

- + B type fully automatic front loading autoclave with a capacity of 45 Ltr / 29Ltr
- + With vacuum, dry cycle & micro processor controlled
- + Touchable screen makes it easy for operation
- + User -friendly operation Interface with elegant design
- + Circuit Breaker for Cuts off external power supply during device malfunction to prevent fire
- + Safety Valve for releases pressure overload and restores idle chamber pressure when the chamber is over-pressured
- + Air Filter for Filters airborne impurities into the chamber to prevent contamination
- + Two Vacuum Pumps for Sufficient vacuum with superior drying results.
- + Two condensers for faster cooling, prolongs life time
- + Two types of record (Dural Recording) - USB & Micro Printer (optional)



Available in 29Ltr also