

XIOS XG Select XIOS XG Supreme



- ✦ CMOS - APS Technology (Breakthrough in Technology)
- ✦ Unique replaceable cable
- ✦ Faster sensor-change
- ✦ Flexible expansion
- ✦ Brilliant HD image quality with 33.3 LP/MM resolution with XG Supreme (Highest in the world)
- ✦ Wi-Fi, cable-free transmission (Optional) for easy room-change with sensor
- ✦ Complete set of positioner

5 Never Before
Years Warranty

CE ISO 9001:2008 ISO 13485:2012 + AC:2012

Made in Germany

1
GLOBAL NUMBER


Simple cable replacement



HD image quality



Flexible expansion



| Specifications | Sensor Size 0 | Sensor Size 1 | Sensor Size 2 |
|--------------------|---|---------------|---------------|
| Technology | CMOS APS | | |
| Greyscale | Select 4096, Supreme 16384 | | |
| Outside dimensions | 32x23 mm | 25.5x38.9 mm | 31.2x43.9 mm |
| Thickness | Select 6.3mm, Supreme 7.5mm | | |
| Active sensor area | 24x18 mm | 30x20 mm | 36x25.6mm |
| Interface | Select - High speed USB 2.0, Supreme - WiFi | | |
| Cable length | Select - 2.7 mtr + 2 mtr, Supreme - 45 cm | | |
| Sensor life | Greater than 4,00,000 exposures | | |