

XIOS XG Select XIOS XG Supreme



- ✦ CMOS - APS Technology (Breakthrough in Technology)
- ✦ Unique replaceable cable
- ✦ Faster sensor-change
- ✦ Flexible expansion
- ✦ Brilliant HD image quality with 33.3 LP/MM resolution with XG Supreme (Highest in the world)
- ✦ Wi-Fi, cable-free transmission (Optional) for easy room-change with sensor
- ✦ Complete set of positioner

5 Never Before
Years Warranty

CE ISO 9001:2008 ISO 13485:2012 + AC:2012


Made in Germany

1
GLOBAL NUMBER

Simple cable replacement



Flexible expansion



HD image quality



Specifications	Sensor Size 0	Sensor Size 1	Sensor Size 2
Technology	CMOS APS		
Greyscale	Select 4096, Supreme 16384		
Outside dimensions	32x23 mm	25.5x38.9 mm	31.2x43.9 mm
Thickness	Select 6.3mm, Supreme 7.5mm		
Active sensor area	24x18 mm	30x20 mm	36x25.6mm
Interface	Select - High speed USB 2.0, Supreme - WiFi		
Cable length	Select - 2.7 mtr + 2 mtr, Supreme - 45 cm		
Sensor life	Greater than 4,00,000 exposures		