

DU32L + DS08 + DLG101



Portable Dental Chair with Unit & Stool



ISO 9001:2008 ISO 13485:2012 + AC:2012

DU 892

PORTABLE DENTAL CHAIR

- ✦ Light weight & sturdy
- ✦ Comfortable & safe for up to 100 kg
- ✦ Back adjusts from upright to spine position
- ✦ Adjustable seat height from 8" to 22"
- ✦ Comes in strong padded case
- ✦ Folded dimensions : 23"X34"X8"
- ✦ Retractable handle for convenient transportation

MOBILE UNIT INCLUDES (BRIEFCASE)

- ✦ 1-Saliva Ejector
- ✦ 1-3 Way Syringe
- ✦ 1-Air Rotor Point
- ✦ 1-Air Motor Point
- ✦ 1-Built in Air Compressor
- ✦ Weight 26 Kg
- ✦ Voltage(V) : 110-240
- ✦ Power(W) : 550
- ✦ Tank(L) : 4
- ✦ Noise Level(dB) : ≤60
- ✦ Dimensions(mm) : 450 X 310 X 680

