

STRONG STRAIGHT HANDPIECE

- ✦ Speed 35,000 RPM
- ✦ Light weight - 56 gm
- ✦ 1:1 ratio handpiece
- ✦ Autoclavable
- ✦ HP Bur \varnothing 2.35 mm

STRONG CONTRA HANDPIECE

- ✦ Speed 30,000 RPM
- ✦ Light weight - 46 gm
- ✦ 1:1 ratio handpiece
- ✦ Autoclavable
- ✦ HP Bur \varnothing 2.35 mm

Handpieces

