

THE LATEST DIGITAL INTRA ORAL SENSOR FOR MODERN IMAGING SOLUTION

THE NEW FRONTIER OF DENTAL IMAGING

Digital dental imaging technology at an unbelievable price. Get ready for crystal clear dental images, which will help you in faster diagnosis and result in happier patients.

Discover a whole new level of clarity.

Key Features :

- ✦ CMOS with APS technology
- ✦ Available sensor sizes - 0, 1 & 2
- ✦ Round edges for patient's comfort
- ✦ Superior image quality
- ✦ User friendly software
- ✦ Essential in successful implant surgery and endodontic treatment
- ✦ Superior image quality with reduced scan time with low X-Ray dose or radiation
- ✦ Plug-in type USB 2.0 PC guarantees the user convenience and simple interaction
- ✦ Flexible cable - Bend resistant and strain relief
- ✦ ORISWIN DG SUITE integrated program for diagnosis, communication and patient database tool with implant simulation



Specifications	Sensor Size 0	Sensor Size 1	Sensor Size 2
Technology	CMOS APS		
Greyscale	4096		
Outside dimensions	31x22 mm	37x24 mm	43x30 mm
Thickness	<5mm		
Active sensor area	24x18 mm	30x20 mm	36x25 mm
Interface	High speed USB 2.0		
Cable length	1.8 mtr + 5 mtr		
Sensor life	Greater than 4,00,000 exposures		