

IMAGE PLATE SCANNER WITH TOUCHSCREEN

The Image Plate Scanner with touchscreen for intraoral formats

The new VistaScan Mini View image plate scanner enables the intuitive, efficient and time-saving digitisation of image plates. Amongst other things, its large touchscreen with its easy-to-use user interface contributes to this. Its compact size and integrated wireless functions make the device really flexible.

This makes image plate diagnostics for dentists quicker, more reliable and more convenient with the VistaScan Mini View.

Digital X-Ray scanner for size 0, 1, 2, 3, 4 image plates. Advanced model with a **touchscreen and WIFI connectivity**.

Easy and cost-effective upgrade to high resolution, digital X-Ray imaging for the whole practice

Features

- ✦ Top image quality
- ✦ High-resolution touchscreen 4.3"
- ✦ Scan Manager for optimum workflow
- ✦ For all intraoral formats (0, 1, 2, 3, 4)
- ✦ Reliability thanks to internal memory
- ✦ PC connection via WLAN/LAN
- ✦ Standalone mode possible
- ✦ Penta Prism technology for scanning

Vistascan Mini View Compact, Fast & Efficient...!!



Made in Germany



Excellent image quality. Thanks to PCS technology

Best in class



reddot design award
winner



NOMINIERT

Data & Facts

VistaScan Mini View	
Screen	4.3" touchscreen, 800 x 480 pixels, 16.7 million colours
Effective resolution (lp/mm, dpi)	22 (1100 dpi)
Theoretical resolution (IP/mm, dpi)	40, 2000
Weight (kg)	Approx. 7 kg
Dimensions (H x W x D mm)	275 x 226 x 243
Standby function	Yes
Interfaces	LAN, wireless