


- ✦ Independent up & down movement and backrest movement of the chair which can be adjusted by the doctor as per convenience
- ✦ Double articulated headrest for handicapped & pediatric patients
- ✦ Right arm rotatable at 90 degree, for easy entry & exit of the patient
- ✦ 7 programs in the chair: 4 programmable working positions, other are spitting, last working position & auto return to zero
- ✦ Chair with anticrushing & movement locking system
- ✦ Integrated bio-system for disinfection of handpiece tubings internally after every treatment
- ✦ Doctor's stool with adjustable backrest tilt & height adjustable, foot ring for attaining a correct working posture by the dentist especially during long dental treatment
- ✦ 3 axis movement of Cold & White 5 LED touchless sensor light (35,000 Lux) with variable intensity
- ✦ Removable & autoclavable handles for dental light head
- ✦ With original monitor arm



CE ISO 13485:2012 + AC:2012 / ISO 9001:2008

MULTIFUNCTIONAL FOOT CONTROL

- ✦ Emergency red switch to lock the chair in working position in case of any accidental obstruction
- ✦ One switch for 4 working position, one switch for last working / spitting position, one switch for auto return to zero, two independent switches for seat & backrest movement
- ✦ One switch for reversal of micromotor, one switch for light on / off
- ✦ Joy stick foot pedal for activation of air rotor handpiece with chip blower / water cut facility

ASSISTANT ARM

- ✦ 2 points - one for low vacuum suction (saliva) and other for high vacuum suction Cannulas only (blood + other debris)
- ✦ 1 point for 3 way syringe

OVERHEAD DELIVERY UNIT

- ✦ With Pneumatic locking arm to prevent any unnecessary displacement of the Unit
- ✦ Feather touch control panel for easy touch play operations
- ✦ 3 points for air rotor with option for fiber optic handpiece
- ✦ An option for inbuilt scaler for oral prophylaxis
- ✦ One 3 way syringe
- ✦ Removable and autoclavable silicone pads for complete germ free operations
- ✦ A pressure gauze on the side of the delivery unit
- ✦ With Inbuilt LED X-Ray viewer

WATER UNIT

- ✦ 180 degree movable water unit
- ✦ Removable spittoon for easy cleaning
- ✦ Dedicated switches for triggering water flow in spittoon & cup filler

