

**The Image Plate Scanner
for all Intraoral Sizes**

The VistaScan Mini image plate scanner makes image plate diagnostics even faster for dentists. The compact device is particularly easy to use and requires a minimum of space – so that it can be installed in the treatment room. The advantage: X-ray and scanning directly at the chairside with full flexibility in the image formats. The reusable VistaScan image plates are read out in top quality within seconds. There has never been a better time to change over to image plates.

Digital X-Ray scanner for size 0, 1, 2, 3, 4 image plates. Easy and cost-effective upgrade to high resolution, digital X-Ray imaging for the whole practice

Features


- ✦ For all intraoral formats
- ✦ Highest image quality
- ✦ Flexible image plates
- ✦ Compact design
- ✦ Ideal chairside appliance
- ✦ Sophisticated operating concept
- ✦ PC interface via USB or network
- ✦ Penta Prism technology for scanning
- ✦ iPad supported App (Optional)

Data & Facts

VistaScan Mini Plus	
Display	Yes
Plate sizes	0 to 4
Effective resolution (lp/mm, dpi)	22 (1100 dpi)
Theoretical resolution (lp/mm, dpi)	40 (2000 dpi)
Weight (kg)	6.5 kg
Dimensions (H x W x D mm)	226 x 234 x 243
Standby function	Yes
Interfaces	USB/Network



Vistascan Mini Plus
compact without compromise



The effective image resolution of X-ray systems in comparison

The VistaScan systems from Dürr Dental achieve the highest resolution compared with the competitor. Latest studies show that image plate technology is the only sure way to uncompromised digitisation