

### X-RAY MACHINE


- ✦ Option to choose between Floor & Wall Mount (AC model)
- ✦ Soft Scissor arm
- ✦ 70 kVp/7mA X-Ray tube
- ✦ Focal Spot 0.8 mm
- ✦ Focal Length 20 cm
- ✦ Inherent Filtration  $\geq 0.5$  mm AL
- ✦ Exposure Time 0.06 to 2.0 second
- ✦ Tube Head & cone are internally LEAD coated
- ✦ Remote control for exposure
- ✦ Microprocessor controlled

#### RAY68(M)

Floor Mount



CE ISO 13485:2003




#### RAY68(S)

Wall Mount

### DC X-RAY MACHINE

- ✦ Constant potential high frequency (DC) generator
- ✦ Very small focal spot (0.4 mm) capable of providing sharp, detailed images
- ✦ Soft Scissor arm
- ✦ Tube Head - Toshiba
- ✦ 65 kVp/7mA X-Ray tube
- ✦ Focal Length 20 cm
- ✦ Total Filtration  $\geq 2.0$  mm AL
- ✦ Exposure Time 0.01 to 2.0 second
- ✦ Tube Head & cone are internally LEAD coated
- ✦ Microprocessor controlled
- ✦ Separate exposure switch



#### Ray 68 DC

Wall Mount

