

**ASPIRE**



- ✦ Independent up & down movement and backrest movement of the chair
- ✦ Proximity sensor based LED light with variable intensity 8,000 - 20,000 Lux
- ✦ DC Motor
- ✦ Programable Dental Chair
- ✦ Doctor's stool with adjustable height
- ✦ Normal balance arm system

**MOVEABLE FOOT CONTROL**


- ✦ Fixed on chair base for up & down
- ✦ Independent foot control for air & water

**ASSISTANT ARM**


- ✦ Movable assistant arm
- ✦ Feather touch panel for all controls of the chair
- ✦ 2 suction point for low & medium suction
- ✦ 1 - 3 way syringe
- ✦ 1 - Extra holder

**DELIVERY UNIT**

- ✦ Pneumatic locking Arm
- ✦ Wide instrument tray
- ✦ Feather touch control panel for all controls of chair
- ✦ 1 point for 3 way syringe
- ✦ 3 points for air rotor
- ✦ A Pressure Gauge present for the air rotor
- ✦ 1 LED X-Ray viewer




#A1  
Sky Blue




#A3  
Light Green



#A6  
Yellow



#A9  
Dark Blue



#A10  
Apple Green



#A11  
Orange



#A12  
Pink

**WATER UNIT**

- ✦ Ergonomic design water unit
- ✦ Removable & easy to clean ceramic spittoon
- ✦ Low & high vacuum suction
- ✦ Cup filler with normal water
- ✦ Programmable water flow in cup filler
- ✦ Built in water purified system