

PEADO DENTAL CHAIR

- ✦ DC Motor
- ✦ Independent Up & Down Movement
- ✦ LED Lamp
- ✦ Right arm movable in up & down position for easy seating
- ✦ TFT screen with memory card for playing videos to attract patient's attention
- ✦ Multifunctional foot control

WATER UNIT

- ✦ Cartoon cuspidor with water canulas
- ✦ Automatic heating-water supply system
- ✦ Option for normal & warm water

CARTOON FISH SHAPE DELIVERY UNIT

- ✦ Big LED X-Ray viewer
- ✦ 1 point for 3 way syringe
- ✦ 2 points for air rotor
- ✦ 1 point for air motor

ASSISTANT UNIT

- ✦ Feather touch control panel
- ✦ 2 Suction point for low & medium vacuum suction
- ✦ 1 Optional point



PEADO+

**Dental Chair for Children
Kids Love it**



EC

**A Special Chair for
your Special Customers "Children"**



Dental Chair

Autoclave

Compressor & Suction

Endomotor & Apex Locator

Implant Micromotor & Piezosurgery Unit

Clinical Equipment