

AUTOCLAVE

BioClave 12L

- ✦ N-type semi automatic Autoclave
- ✦ Front loading bioclave with a capacity of 12 Ltr
- ✦ Microprocessor controlled
- ✦ Ensures biosafety & Sterilization of non porous instruments
- ✦ Completes dry cycle with semi - opened door
- ✦ Door lock for safety
- ✦ Stainless steel chamber and 2 stainless steel trays for surgical instruments
- ✦ Silicone ring seals the chamber
- ✦ Thermostat prevents over heating
- ✦ Diameter of the chamber 21.4 cm
- ✦ Dimension - 40.7 cm width x 38.5 high x 46.5 cm deep
- ✦ Net weight 23 Kg

AUTOCLAVE VITALE 12 N TYPE

- ✦ Digital pressure gauge with LED bar
- ✦ Unique Program of Sterilization
- ✦ Capacity 12 Ltr
- ✦ Keyboard control in turquoise
- ✦ Chamber with 2 trays in anodized aluminum
- ✦ Drying with door closed jar
- ✦ 16 security systems
- ✦ System microcontroller
- ✦ Output Power 1600W
- ✦ Chamber dimension : (W x D) 220 x 330mm
- ✦ Autoclave dimension (W x L x H) : 335 x 335 x 480mm
- ✦ Net weight 24 kg

QI 15L / SEA 15L-N-LED

- ✦ N type fully automatic front loading autoclave with a capacity of 15 Ltr
- ✦ With dry cycle autoclave & micro processor controlled
- ✦ Wrapped / unwrapped, non porous, solid instrument can be sterilized
- ✦ Digital display of pressure and temperature
- ✦ Inbuilt Sensor for distilled water
- ✦ Double safety lock on the door
- ✦ Silicon Gasket prevents any leakage of air
- ✦ Diagnostic automatic alarm for any malfunctioning
- ✦ Can be used at 121°C / 134°C during sterilization
- ✦ 2 Storage tanks for fresh and used water
- ✦ Chamber size 230 mm x 360 mm
- ✦ Net weight 50 Kg
- ✦ No. of Tray 3 Nos.

GNATUS



ISO 13485:2012 + AC:2012 / ISO 9001:2008

Cristófoli



CE EN ISO 13485:2012 + AC:2012

Runyes



15 Ltr N-Type front loading Fully automatic

CE ISO 13485:2003