

FENG 22L / SEA 22L-B-LED

- ✦ B type fully automatic front loading autoclave with a capacity of 22 Ltr
- ✦ With vacuum, dry cycle & micro processor controlled
- ✦ Wrapped / unwrapped, porous, solid instrument including cotton swabs can be sterilized
- ✦ Bow & Dick test qualified, Inbuilt Sensor for distilled water
- ✦ Digital display of pressure & temperature
- ✦ Double safety lock on the door
- ✦ Silicon Gasket prevents any leakage of air
- ✦ Diagnostic automatic alarm for any malfunctioning
- ✦ Can be used at 121°C / 134°C during sterilization
- ✦ 2 Storage tanks for fresh & used water
- ✦ Chamber size 250 mm x 450 mm
- ✦ Net weight 55 Kg
- ✦ No. of Tray 3 Nos.



22 Ltr B-Type front loading Fully automatic

CE ISO 13485:2003

SEA-22L-B-LED Touch Autoclave

- ✦ B class European standard, 3 times prevacuum autoclave with a capacity of 22L
- ✦ MOTOROLA sensor
- ✦ Italian water pump
- ✦ Intelligent Automatic control System
- ✦ Automatic electric door system
- ✦ Automatic cleaning reminder every 100 cycles
- ✦ Open water tank - easy to be cleaned
- ✦ Big capacity 3.5L clean water tank & 5L used water tank
- ✦ Chamber size ϕ 249 x 450
- ✦ Product Size 594mm x 468mm x 453mm
- ✦ Installation Size 60 cm
- ✦ Weight 50kg

NEW



22 Ltr B-Type Touch automatic

CE ISO 13485:2003