

- ✦ Medical grade air
- ✦ Super silent
- ✦ Oil free engine
- ✦ Air outlet filter
- ✦ Thermal protector against overheating



BioQualyAir 32L



BioQualyAir 40L



BioQualyAir 65L

Model	Power	L/min	dB(A)	Bar	Ltr.	Dimensions (mm) (L X W X H)
BioQualy Air 32 L	550W	122	55	6-8	32	-
BioQualy Air 40 L	1100W	212	65	6-8	38	400 x 695 x 400
BioQualy Air 65 L	2200W	424	71	6-8	65	510 x 765 x 510

ISO 13485:2012 + AC:2012 / ISO 9001:2008

## BioVacSec Motorized Suction

### NO WATER. GREAT PERFORMANCE

**Dry running system. Does not consume water in its operation. Can be installed any where, without the need of water point nearby. Low power consumption**

- ✦ Very low noise level
- ✦ Does not bother the patient nor the dentist
- ✦ It can be installed up to 7 meters from the point of suction, while maintaining its performance
- ✦ Drains by gravity - can be connected directly to the common sewer
- ✦ Power consumption less than 0.45kwh
- ✦ Anti Vibration System with cushion for fixing the engine & rubber feet
- ✦ Easy access to the debris filter separator
- ✦ ½ HP/ 0.37 kw motorized suction
- ✦ Level sensors - Controls the level of the reservoir, avoiding overflow
- ✦ Reservoir with 1 liter capacity & Serves up to 2 dental clinics



ISO 13485:2012 + AC:2012 / ISO 9001:2008