

EMS-200

Four-in-One (Endo Motor, Apex Locator, Filler & Obturator System). Complete Solution for an Endodontist

Features

- ✦ One touch color wide screen
- ✦ Compact & light weight
- ✦ Endo motor & Apex Locator at the same time
- ✦ Options to Mount on dental chair or movable cart type
- ✦ Best accuracy & systematic diagnosis

Endo Motor Features

- ✦ Auto-reverse mode
- ✦ Torque control from 0.6 to 6.4Ncm
- ✦ Long shelf life with swiss motor
- ✦ 10 Memory saving function
- ✦ Reciprocating mode

Obturator (Fill) Features

- ✦ Easy install for fill needle
- ✦ Fixed cartridge needle
- ✦ Fully motorized obturation for smooth filling
- ✦ Long shelf-life with swiss motor
- ✦ Comfortable grip with ergonomic design
- ✦ Auto reverse function when cartridge is empty
- ✦ Automatic off after 5 minutes
- ✦ Option to reused cartridge

Apex Locator Features

- ✦ Best accuracy in any root canal
- ✦ Self test for system failure
- ✦ Simple, easy to use
- ✦ Color display

Obturator (Pack) Features

- ✦ 360° rotating for pack tip
- ✦ Outstanding heating performance (Max. 400°C)
- ✦ Comfortable grip
- ✦ Hold switch for safety
- ✦ Temperature range (140~400°C) displayed in LED window
- ✦ Automatic off after 5 minutes

A DREAM MACHINE
for the real Endodontist in you.



Dental Unit For Endodontic Treatment

4 IN 1 SYSTEM

ENDO MOTOR

Shape the Root Canal

APEX LOCATOR

Measure The Working Length

OBTURATOR PACK

Pack Root Canal

OBTURATOR FILL

Fill Root Canal