

ENDO-MATE DT

ENDO MICROMOTOR WITH TORQUE CONTROL AND AUTO REVERSE

- ✦ Variable Speed
- ✦ 9 programmable memory
- ✦ 6 head angles for ease of treatment
- ✦ Push type Ultra miniature head
- ✦ On/Off switch on hand piece
- ✦ Wide screen provides high visibility
- ✦ Compact desktop unit for all types of rotary files
- ✦ Speed range : 100-13,000 min⁻¹
(With 20:1, 4:1, 1:1, 16:1 geared heads)
- ✦ Dimensions : W92 x D148 x H124 (mm)
- ✦ Weight : Unit 456g/Motor handpiece & cord 92g
- ✦ Torque Level : 0.8-6.5 Ncm (When using 20:1 Handpiece)
- ✦ Battery / AC operated

Endodontics has never been so EASY...!!



ENDO-MATE AT

STYLISH & COMPACT ENDO MICROMOTOR

- ✦ Fits all major brands of Ni-Ti Rotary Files
- ✦ Torque control, with integrated auto-reverse function
- ✦ Compact Stylish Control/ Monitor pad
- ✦ Large screen Size for immediate visibility
- ✦ Flag control panel for simple & user friendly operation
- ✦ Light weight, Smart hand piece offering excellent visibility
- ✦ Convenient On/Off Switch
- ✦ Power supply AC 230 V 50/60 Hz
- ✦ Nos of programs 9
- ✦ Max. torque 7 Ncm (with 20:1 head)
- ✦ Speed range 100-13,000 min⁻¹
- ✦ Dimensions W268 x D230 x H130 (mm)
- ✦ Weight 342 gms only (Including motor handpiece & cord)
- ✦ AC operated



ENDO-MATE TC2

CORDLESS ENDODONTIC HANDPIECE WITH TORQUE CONTROL & AUTO REVERSE

- ✦ Simple 5-key operation
- ✦ Light weight
- ✦ Cordless unit
- ✦ Alarm function alerts
- ✦ Variable torque
- ✦ Variable speed
- ✦ 5 programs
- ✦ Max. torque : 4 Ncm (when using 20:1 head)
- ✦ Speed range: 100-2,250 min⁻¹ (With 20:1, 16:1, 10:1, 4:1 geared heads)
- ✦ Charge/continuous use : 90 min/1.2h (depending on workload)

