

MPAS HEAD

You can use this with Endo-Mate and Apex Locator simultaneously to take benefit in **RCT**



iPexII

**PRECISION DIGITAL APEX LOCATOR
FEATURING THE NEXT-GENERATION
SMARTLOGIC CONTROLLER**

- ✦ High-precision root canal measurement generated by using SmartLogic
- ✦ Auto-detects the apex accurately in any conditions, dry or wet
- ✦ A stylish and compact body requiring only a small space. Can be located on the dental equipment console
- ✦ Three different alert sounds depending on the location of the file tip, ensuring precise procedures
- ✦ A big 3 colour LCD panel with very low reflectivity, allows Clinicians to accurately & instantly monitor procedures
- ✦ Energy saving & low operation cost
- ✦ Great accuracy in wet or dry canal
- ✦ The NSK iPex has multi-frequency measuring technology
- ✦ User friendly interface
- ✦ Large Hi-contrast LCD screen
- ✦ Audible warning system
- ✦ Auto power shut-off
- ✦ Battery : Continuous operation for 60 hours
- ✦ Display : 3 color LED indicating of file movement inside the canal
- ✦ Power supply : for 1.5V AAA (LR03 or R03) batteries
- ✦ Continuous use : Approx. 60 hours
- ✦ Dimensions : W60 x D60 x H86 (mm)
- ✦ Weight : 76g (unit without batteries)

Compact and Stylish

Lightweight and with a base just 6 cm in diameter, the iPexII can quickly be set up and used in virtually any environment. The LCD panel is pitched at the ideal angle that allows Clinicians to easily & without distraction monitor conditions ensuring safe and smooth root canal procedures

The display panel is easily viewable from an operating stance & also minimizes reflectivity, while maintaining the Clinicians focus on the procedure.



Dental Chair

Autoclave

Compressor & Suction

Endomotor & Apex Locator

Implant Micromotor & Piezosurgery Unit

Clinical Equipment

CE ISO 9001:2008 ISO 13485:2003