



The next instrument for your tray
The Laser for every Dentist, every day

2 Years Warranty

Soft Tissue Diode Laser 940nm-5w



FDA ISO 13485:2012 + AC:2012

With no foot pedal, power cord or external controls,
the revolutionary iLase™ takes dental lasers beyond portable.

Main Features

- + Finger switch Laser Actuation
- + Maximum of 5W Peak output Power & Continuous Wave Output Power
- + Factory loaded pre-set values
- + Battery Power
- + 940nm Wavelength
- + Comfort Pulse™
- + Bendable, disposable, multiple length, multiple diameter tips
- + User Friendly display & navigation

**Powerful, Affordable,
Beyond Portable**

iLase™ Technical Specifications

Handpiece	
Dimensions	Length : 7.2" (183mm), with battery attached Diameter : 0.74" (18.7 mm)
Weight	0.22 lb (98g), with battery attached
Max. Peak Output Power	5 W
Wavelength	940 ± 15nm
Three Power Modes	Continuous Wave ComfortPulse™ 1 ComfortPulse™ 2
Battery Power	Single, rechargeable Li-Ion assembly; 3.7 VDC, 650 mA-h
Presets	10 factory-loaded & user-customizable 2 extra user-customizable
CHARGING STATION	
Dimensions	W x H x D (4.72" x 2.5" x 3.95") (120 x 70 x 100 mm)
Weight	0.90 lb (412g)
Charging Ports	4