

VarioSurg3

Optic

Exceptional Power for Faster, More Accurate Procedures

The exceptional power VarioSurg3 that enables precise, stressfree surgical procedures even when approaching hard bone in delicate areas.

The compact body offers a wide variety of functions. The large backlit LCD display & intuitive controls help reduce operating errors & make procedures safer.

Applications


- ✦ Bone surgery
- ✦ Sinus lifts
- ✦ Extraction
- ✦ Implant preparation

Main Features

- ✦ Major power boost
- ✦ 3 Modes (SURG, ENDO, PERIO)
- ✦ Selectable burst mode
- ✦ 9 Programs (SURG x 5, ENDO x 2, PERIO x 2)
- ✦ Adjustable LED Light Intensity
- ✦ Automatic cleaning mode
- ✦ Feedback function
- ✦ Auto tuning function
- ✦ Dynamic link feature major power boost
- ✦ LED handpiece with 2 meter cord
- ✦ Frequency : 28 – 32 kHz
- ✦ Irrigation flow rate : 10 – 80 mL/min
- ✦ Dimension : W265 X D220 X H100 (mm)
- ✦ Power supply : AC120V or 230V 50/60Hz
- ✦ Maximum output 25 W
- ✦ Multifunctional foot control (FC-78)
- ✦ Sterilization cassette
- ✦ Basic H-S - Kit : (H -SG1, SG3, SG5, SG6D, SG7D, SG11 & Tip Holder)
- ✦ Handpiece stand, Irrigation tube (5 pcs) & other accessories



50% Increased Power Focus on the Procedure



NSK's twin LEDs deliver natural daylight quality illumination & significant advantages to the way you work

Clearer : LEDs generate natural daylight-quality light, which ensures clearer vision than halogen light sources.

More Durable : LED light sources are longer lasting & more durable than halogen bulbs.