




PIEZO-SURGERY UNIT

The Piezo Surgery Unit from Woodpecker combines control with ease of use in one handy package. And with an LED handpiece to ensure better vision, you have absolute precision. Use this unit to perform bone cutting with more control and precision, facilitating tooth extraction and bone surgeries.

- ✦ Easy operation: Separate functions for bone root and cleaning
- ✦ Enhanced operator control and sensitivity leading to better feel and precision for the cutting action
- ✦ Safe cuts
- ✦ 14 Tips to various procedure
- ✦ Comes with 2 Handpieces
- ✦ Removable peristaltic pump
- ✦ Power Presets: According to Bone density
- ✦ Better irrigation flow control



Bone Mode - Voltage signal with high and low amplitude alternation adjusted period supplies high-power energy for the handpiece, ensuring high cutting efficiency even when cutting the hardest texture.


Perio&Endo Mode - Sine voltage signal with constant amplitude, supplies the handpiece with stable power, ensuring precise treatment.

5 Years Warranty on PCB

- Now Available with
- ✦ +60 Watts
 - ✦ +7 Additional Tips
 - ✦ +One Extra handpiece with cord

Working 10 min under full power, test the temperature of handpieces

The temperature below is the maximum temperature



Well-known European Brand S 30°C

DTE Brand 33°C

Well-known European Brand M 33°C

Well-known European Brand D 38°C

50°C 26°C Color temperature value

Measurement tool : FLIR -E60 infrared camera

