

- ✦ Connect to any standard air tubing coupling with any dental chair
- ✦ No need of circuit board, tubings & bulbs
- ✦ High intensity & stable self LED illumination through mini power generator
- ✦ Cellular glass optic - pure daylight quality illumination
- ✦ Stainless steel : durable, stable & light weight
- ✦ Available with Miniature Head M500LG (380,000 - 450,000 RPM) Head Size: $\phi 10.6 \times H12.4\text{mm}$, 25,000 LUX & Standard Head M600LG (330,000 - 430,000 RPM) Head Size: $\phi 12.1 \times H13.3\text{mm}$, 25,000 LUX
- ✦ Clean head system
- ✦ Push button chuck
- ✦ Quattro Spray

Visibility & Access :

Improved visibility and direct access to the treatment field allows faster and more accurate preparations.

LED Durability :

The LED light source is significantly more durable than quartz. Replicated daylight illumination quality.



Optic - NSK LED



Non Optic


S-Max pico

The World's Smallest Head

The innovative air turbine with ultra mini head and super slim body shape

NSK has developed special ultra short shank burs to operate specifically with the pico handpiece. These high quality precision burs are concentric and durable. Use of NSK ultra short shank burs maximizes the mini advantage of the pico handpiece.

- ✦ Stainless steel body
- ✦ ISB turbine
- ✦ Cellular glass optics
- ✦ Ceramic bearings
- ✦ Clean head system
- ✦ Single spray
- ✦ Power 9 W
- ✦ Speed 360,000~450,000 min⁻¹
- ✦ Push Button chuck



Connect to NSK Couplings

