

Pana·Max PLUS

Non-Optic


- ✦ Stainless Steel Body
- ✦ Ceramic Bearings
- ✦ ISB Turbine
- ✦ Clean Head System
- ✦ Push Button Chuck
- ✦ Quattro Spray
- ✦ Miniature Head Speed : 380,000-450,000 min⁻¹
- ✦ Standard Head Speed : 340,000-410,000 min⁻¹

Ti-Max X

eXtreme Premium Series

- ✦ Titanium Body with Scratch Resistant DURACOAT
- ✦ Cellular Glass Optics
- ✦ Clean Head System
- ✦ Push Button Chuck (X95L/X95, X25L/X25, X15L/X15)
- ✦ Micro Filter (X95L/X95, X25L/X25, X15L/X15)
- ✦ Quattro Spray
- ✦ For FG burs (ø1.6)
- ✦ Ceramic Bearings
- ✦ Max. Speed: 200,000 min⁻¹

- ✦ Titanium Body with Scratch Resistant DURACOAT
- ✦ Cellular Glass Optics
- ✦ Clean Head System (X85L/X85, X10L/X10, X20L/X20)
- ✦ Push Button Chuck (X85L/X85, X10L/X10, X20L/X20, X35L/X35, X70L/X70, X75L/X75, Ti-ENDO)
- ✦ 64:1 Reduction
- ✦ Single Spray
- ✦ For Ni-Ti File
- ✦ Max. Speed : 300 min⁻¹



CE ISO 9001:2008 ISO 13485:2003