

**D7 LED SCALER**

- ✦ Detachable handpiece with light, better visibility, autoclavable
- ✦ Works on 3 modes : Scaling, Perio & Endo
- ✦ Clinical solutions are applicable for automatic water supply mode, improving the performance of perio and endo function dramatically
- ✦ Less heat from tip during operation. Improve comfort degree of patients
- ✦ Output power 3W~20W
- ✦ Frequency 28kHz±3kHz
- ✦ Dimensions 237mm x 206mm x 95mm



Fiber optic (LED) Handpiece

**DTE - D5**

**Optic / Non-Optic**

- ✦ Detachable & autoclavable handpiece
- ✦ Works on 3 modes : General, Perio & Endo
- ✦ Automatic frequency tracking
- ✦ With 5 Titanium Scaling tips (3 S, 1 P, 1 Stain Removal) & 1 Endo attachment
- ✦ Power input : 100V to 240V~ 50Hz/60Hz 1.2A (Max)
- ✦ Output Power : 3W to 20W
- ✦ Frequency : 28 kHz + 3 kHz
- ✦ Dimension : 188 mm X 159 mm X 81 mm
- ✦ Main unit weight : 0.65 kg



Fiber optic (LED) Handpiece

**DTE - D3**

**Optic / Non-Optic**

- ✦ Detachable & autoclavable handpiece
- ✦ Works on 2 modes : General & Perio
- ✦ Automatic frequency tracking
- ✦ With 5 Tips (3 S, 1 P, 1 Stain Removal)
- ✦ Power input : 100V to 240V~ 50Hz/60Hz 1.2A (Max)
- ✦ Output power : 3W to 20W
- ✦ Frequency : 28kHz + 3kHz
- ✦ Dimension : 196 mm X 134 mm X 80 mm
- ✦ Main unit weight : 0.62 kg



Fiber optic (LED) Handpiece