

INBUILT SCALER V2

- ✦ Detachable handpiece
- ✦ 2 modes
- ✦ With 5 Tips (3 S, 2 P)
- ✦ With Torque wrench
- ✦ Power input : 100V to 240V~ 50Hz/60Hz 1.2A(Max)
- ✦ Output Power : 3W to 20W
- ✦ Frequency : 28kHz + 3kHz
- ✦ Dimension : 74 mm X 56 mm X 38 mm
- ✦ Main unit weight : 0.2 kg
- ✦ Inbuilt Scaler V1 LED & V3 LED also available



CE ISO 13485:2012 / AC:2012

LIGHT CURE UNIT

- ✦ Strong optical output power
- ✦ Three treatment modes : Continuous, Ramp & Pulse
- ✦ Digital time setting : 0~20s for LED curing light
- ✦ Full color LCD display
- ✦ LED cold light output can treat all dental tissue
- ✦ Alarm for low voltage

Technical Specifications

- ✦ Working voltage : 100-220 V AC, 50-60 Hz
- ✦ Working volt : 4.2-4.5 volts
- ✦ Electric current : 0.7 amperes
- ✦ Luminous intensity : > 1200 mw/cm² 5W
- ✦ Wavelength : 470 nm
- ✦ Weight : 160g
- ✦ With 1 black probe



CE ISO 13485:2012 / AC:2012

DIGITAL LED

LIGHT CURE UNIT

- ✦ Material of shell handle : Aluminum
- ✦ Size : Diameter 23×160mm
- ✦ Weight : 160g
- ✦ Grade : Ordinary B Type Intermittent working
- ✦ Working Voltage : 100-220V AC, 50-60HZ
- ✦ Grade : II IP40
- ✦ Mode : Three kinds of working mode (Alterability)
- ✦ Working Volt : 4.2~4.5 volts
- ✦ Electric current : 0.7 amperes
- ✦ Battery Type : Lithium Ion
- ✦ Battery capacity : 1400Mah
- ✦ Wave Length : 420-480nm
- Luminous Intensity : 1200mW/cm² 5W

(According to IEC60601-1, curing lights belongs to Class II, I P40)



TULIP