

INTRA ORAL CAMERA

INTRA ORAL CAMERA WITH MONITOR

- ✦ Camera features a 1/4" High Definition image sensor
- ✦ 17" Monitor features a High Definition resolution of 1280 X 1024
- ✦ 8 GB Memory Card (included) allows for easy image transfer & file management
- ✦ Camera buttons are easily accessible to control capture & playback functions
- ✦ Camera is illuminated with six bright LEDs, giving the best contrast, definition & resolution quality
- ✦ Playback format includes : MP3 audio, MP4 video & JPG images
- ✦ Adjustments on monitor include : Contrast, brightness and color
- ✦ Monitor comes equipped with a medical grade tempered glass front panel to minimize scratches
- ✦ Easy integration of LCD monitor with camera on the dental delivery unit due to no additional cable needed
- ✦ Input supports USB, VGA, HDMI & AV



ADVANCE CAM IOC VGA + USB + AV MODEL

- ✦ Intra Oral camera with 6 LED
- ✦ CCD Sensor
- ✦ Memory - Camera 50 pic & USB 10,000 pic
- ✦ VGA & USB 2.0 signal Output
- ✦ Resolution : 2 mega pixels
- ✦ Compatible to USB, TFT Monitor without CPU
- ✦ Lens 0.02 inch pin lens
- ✦ Focus Range 8-12 mm
- ✦ Backup - USB device backup
- ✦ Display modes - Full screen & 4 split modes
- ✦ Intelligent function - Prevent dithering, Anti-alias edge
- ✦ Now available with Window 7 (Compatible)

