

## INTRA ORAL CAMERA

### DIRECT USB IOC


- ✦ Resolution : 2 mega pixels, perfect image
- ✦ Optical Format : 1/4-inch (4:3)
- ✦ Interface : USB output
- ✦ Dynamic Range : 60dp
- ✦ Resposivity : LUX-see (550nm)
- ✦ Plug & play - work immediately
- ✦ 2 blue LED doe carries detection

### WIFI INTRA ORAL CAMERA

- ✦ End-user : Support 4 users at the same time
- ✦ Available Distance :  $\geq 30m$ (open space)
- ✦ Battery capacity : 2100mA/h
- ✦ Working time : 140mins
- ✦ Pixel : 1 million
- ✦ Light Resource : 6pcs white LED light(5600K)
- ✦ Dimming : Five brightness adjustment
- ✦ Focus Range : 10mm-50mm
- ✦ Photography Angle :  $105^\circ$
- ✦ Connection : WiFi
- ✦ Materials : Imported medical resin
- ✦ Main Function : Real-time watching, oral examination
- ✦ Two-way take photo and video
- ✦ Compatible with iphone & Android devices

### X-RAY FILM VIEWER

- ✦ LED light source, working life over 100000 hours
- ✦ Can mount on the wall or place on the desk with 3 positions
- ✦ Fashionable design
- ✦ Intensity of light :  $>3000$  LUX
- ✦ Dimensions : X-Ray viewer - 35,7 x 24,3 x 3,5 cm (W x H x D) / Screen - 30 x 20 cm (W x H)
- ✦ Illuminance : 7.000 Lux
- ✦ Light : Cold, polarized, no oscilation
- ✦ Power : 10 Watts
- ✦ Voltage : 220 V - Automatic voltage selection
- ✦ Net Weight : 1,1 kg



HD : 720P

Cristófoli



Luna

EN ISO 13485:2012 + AC:2012