

## DIGITAL WHITE STAR BLEACHING UNIT

- + Easy to adjust at any angle position
- + 6 pieces high power LED with fan cooling
- + 1-30 min adjustable timer
- + Wide LED display screen
- + Wavelength : 420 nm - 490 nm
- + Light intensity : 3 intensities  
Max>6000 mw/cm<sup>2</sup>
- + Input power : 100-240V AC, 50/60Hz
- + Available in floor & chair mount options



Chair Mount

ISO 13485:2012 + AC:2012

## LED BLEACHING UNIT

- + Flexible arm gives user ease of use
- + Works well on all types of stains - internal stains, tetracycline stains, flourosis
- + Excellent results in single visit
- + Wavelength : 450 - 480 nm
- + Intensity : 8000 - 10000mw/cm<sup>2</sup>
- + 8 pc high LED
- + Available in floor & chair mount options



Floor Mount

## U V CHAMBER

- + S.S. Chamber with 12 S.S Trays
- + Imported UV tube
- + Emits UV light
- + Effective tool for avoiding contamination of instruments once sterilized or autoclaved
- + Can maintain & preserve instruments for days without any fear of contamination
- + Door switch interlock prevents emission of UV light as soon as the chamber is opened to avoid any harmful effect on operator
- + Electronic circuit for instant start of tube

### Technical Specification

- + Chamber Dimension : A - 12 Trays - 300 X 190 X 510mm  
B - 15 Trays - 300 X 190 X 572mm
- + Power Consumption : 8 & 16watts UV Tubes.
- + Power Supply : 220V 50Hz

**NEW**

