

## ULTRASONIC CLEANER

- ✦ Tank Material : SUS304 Stainless Steel
- ✦ Tank Capacity : 10 Ltr
- ✦ Time setting : 1-60min adjustable
- ✦ Temperature setting : 200° C - 800° C adjustable
- ✦ High-performance transducer
- ✦ CE/FCC/RoHS tested and approval
- ✦ More quiet to operate in the hospital
- ✦ Special groove design for better waterproof



CC GB/T/9001-2008 ISO 9001:2008

## ULTRASONIC CLEANER

- ✦ 3-colour LED display : 1-30 min digital timer, 5 temperature settings, working status
- ✦ Advanced programming : Temperature display, Solution usage timer
- ✦ Tank capacity : 3 Ltr
- ✦ Overheat protector : To prevent overheating & life of heater
- ✦ Dual heat protections : Heater turns off after 90 mins for safety
- ✦ Degas features : Removes bubbles quickly
- ✦ Tank size : 25 cm x 17.4 cm x 7.8 cm
- ✦ Power : 150W (AC 220 ~240V 50 Hz)
- ✦ Digital timer : 1-30 mins full range timer
- ✦ Ultrasonic frequency : 36,000Hz
- ✦ Weight : 3.6 kgs



## DIGITAL ULTRASONIC CLEANER

- ✦ 3 color LED display
- ✦ 1 to 30- minute full range timer
- ✦ 5 temperature settings
- ✦ Solution usage timer
- ✦ Degas feature
- ✦ Independent power switch and detachable power cord, convenient to use
- ✦ Tank capacity 3 Ltr
- ✦ Frequency 35,000 Hz
- ✦ Tank Size 25 x 17.4 x 7.8 cm and Weight 5.1 Kg
- ✦ Tank material stainless steel SUS304



### Standard Accessories :



Basket

### Optional Accessories :

- ✦ SS Tray ✦ Rack
- ✦ Beaker ✦ Beaker Holder