

DENTAL SEALER

- ✦ The effective heating control and high quality performance because of the automatic temperature adjustment and stand-by system
- ✦ Reminder system prevents the overheating during sealing
- ✦ The smart operation pad allows to choose the exact bag length to save the paper
- ✦ 12mm sealing seam can definitely guarantee the legal requirements of sterilization
- ✦ Rated power 80W, creating the new record of energy saving in the sealers' family



Galilean, Dental Surgical Loupe (MODEL 2.5 X & MODEL 3.5 X)

- ✦ 2.5 or 3.5 times magnification
- ✦ Offers widest field of view with adjustable focal length
- ✦ Flip up facility
- ✦ Convenient fastener makes it useful for everyone
- ✦ Optimal magnification
- ✦ Also available model GL4-3.5X & GL3-2.5X



Head - Light
for Loupe

AQUA - WATER DISTILLER

Water distiller is the best in its class. This distiller completes one liter of pure steam distilled water in 1 hour. The boil chamber and condensing coils are high grade 304 Stainless Steel. It comes with a high quality plastic collection bottle, the lid and handle are made of food-grade polypropylene, and the lid is removable for convenience. Distilling your own water can't get any easier, simply fill the boil chamber, press the start button and walk away, the distiller will automatically shut-off at the end of the cycle. It's very easy and convenient.

Features

- ✦ Easy to use, produces pure steam distilled water
- ✦ Lightweight compact design, great for home and for medical use
- ✦ Eliminates plastic bottles from the store & saves you money
- ✦ Capacity - 4 Ltr
- ✦ High quality collection bottle for the purest water

