

- ✦ Medical grade air
- ✦ Super silent
- ✦ Oil free engine
- ✦ Air outlet filter
- ✦ Thermal protector against overheating

### BioQualyAir 32L



### BioQualyAir 40L



Model	Power	L/min	dB(A)	Bar	Ltr.	Dimensions (mm) (L X W X H)
BioQualy Air 32 L	550W	122	55	6-8	32	-
BioQualy Air 40 L	1100W	212	65	6-8	38	400 x 695 x 400
BioQualy Air 65 L	2200W	424	71	6-8	65	510 x 765 x 510



### BioQualyAir 65L