

- ✦ Seamless & soft upholstery for patient comfort
- ✦ Synchronized movement of seat & Backrest
- ✦ Adopt ABS engineering plastics and the whole mold making
- ✦ Double articulated headrest for handicapped & pediatric patients
- ✦ Right arm rotatable for easy entry & exist of the patient
- ✦ Touch free (sensor) operation of LED Light with minimum intensity from 8,000- 25,000 Lux
- ✦ 90-degree moveable water unit for 4-hand dentistry
- ✦ 3 customized programs with spitting, last working position & auto return to zero position
- ✦ One-piece PU base, perfect combination of luxury and steadiness
- ✦ All connections in the base of chair
- ✦ Durable & inbuilt foot paddle made of PU material with feather touch control panel on doctor side & assistant side

INNOVA GOLD



Stainless Steel
Foot Control

CE ISO 13485:2003

*Optional



- ✦ When the handpiece is working, dental chair will automatically keep locked to ensure the operation safety
- ✦ Safety with protection against any obstruction when chair goes down
- ✦ Stainless steel base doctor stool
- ✦ Stainless steel progressive foot control

DELIVERY UNIT


- ✦ Three points for air rotor, one 3 way syringe and one optional point
- ✦ Silicon tubing with non-retraction valve
- ✦ Wide instrument tray
- ✦ Fiber optic inbuilt scaler
- ✦ Silicon pad on instrument tray to prevent from corrosion and cross contamination
- ✦ Flip type film viewer

WATER UNIT

- ✦ Water unit have on/off water change, water control buttons
- ✦ Inbuilt water reservoir with empty sensor
- ✦ Temperature control for tumbler water for patient comfort



#C01




#C02



#C03



#C04



#C05

ASSISTANT ARM

- ✦ Rotatable assistant arm with 3 way syringe & low /high suction cannula
- ✦ Detachable filter for easy cleaning
- ✦ Strong high & low vacuum based suction system with time delay function to avoid cross contamination
- ✦ Suction cannula with aspiration control switch.

DCI LIGHT

- ✦ 3 axis rotation (For easy positioning and projection)
- ✦ 3-intensity : 16000 Lux, 26000 Lux & 38000 Lux
- ✦ Precise illumination Pattern : 3" X 6" (Prevents Patient Glare)
- ✦ Color Temperature : 5000K (Accurate shade matching)
- ✦ Color Rendering Index (CRI) : 90
- ✦ Shadow free light (Less eye fatigue & Increase Visibility)
- ✦ Reflective technology : Uses 500 reflective facet (Uniform light pattern)
- ✦ Durable Design : Heavy-duty structure & handle design
- ✦ Hermetically sealed front screen makes external cleaning easier and stops any dust or spring fluid getting inside the light

INNOVA
PAD



Stainless Steel
Foot Control

CE ISO 13485:2003

*Optional



DCI LIGHT

