

### BioVacSec Motorized Suction

#### NO WATER. GREAT PERFORMANCE

**Dry running system. Does not consume water in its operation. Can be installed anywhere, without the need of water point nearby. Low power consumption**

- ✦ Very low noise level
- ✦ Does not bother the patient nor the dentist
- ✦ It can be installed up to 7 meters from the point of suction, while maintaining its performance
- ✦ Drains by gravity - can be connected directly to the common sewer
- ✦ Power consumption less than 0.45kwh
- ✦ Anti Vibration System with cushion for fixing the engine & rubber feet
- ✦ Easy access to the debris filter separator
- ✦ ½ HP/ 0.37 kw motorized suction
- ✦ Level sensors - Controls the level of the reservoir, avoiding overflow
- ✦ Reservoir with 1 liter capacity & Serves up to 2 dental clinics

### MOTORIZED SUCTION

