

DENTAL CHAIR UNIT



**White & Cool
LED Light**



**Modular
Delivery Unit**



Water Unit



**Multifunction
foot control**

DENTAL CHAIR - SPECIFICATION

DENTAL CHAIR

Fully Electrically operated dental chair equipped with state-of-art microprocessor based chair control system

- Feather touch switches for chair operation with Zero programme
- Articulated Headrest with comfortable cushioned rexine
- Swivelling arm rest for easy patient access

DENTAL LIGHT

White & Cool LED Operating light having multi surface reflector with capacity of 30000Lux

- On/Off and intensity control by non-touch sensor with 3 directional movements
- Light Fitted on well balanced arm for operator's convenience.
- Offers an Operating light with removable handle

WATER UNIT

- Porcelain spittoon with auto water connection for cup filler & flush function
- High & Low vacuum motorised suction with auto drain & auto flush system

- Assistant side chair console with High & Low vacuum suction cannulas
- Assistant 3-way syringe with autoclavable tip
- Fast curing Light cure unit

DELIVERY UNIT

- Modular delivery hanging system
- Two air rotor points with water control option
- Micro motor 35,000rpm with digital display of speed
- Non-optic scaler with five scaling tips
- Doctor 3-way syringe
- SS Instrument Tray
- X-Ray viewer

Multifunction foot control can be use for Air Rotor & Micromotor Operating

- Joystick function for chair movement

Dental operator's stool moving on five castors with chromium plated base