

BIOLASE

waterlase^{*}MD TURBO™

An Ideal **All-tissue Laser** for the Dental Specialist

Two Contra-angle 360° Rotatable Handpieces For Maximum Versatility



Extremely Light And Flexible Titanium Fiber Cable Helps Eliminate User Fatigue



Intuitive Operator



A periodontal revolution

- + Deep Pocket Therapy using the patented Radial Firing Perio Tip

Endolase Root Canal Therapy

- + Access cleaning and enlargement and disinfection of the root canal.

Easy Day-to-Day Cavity Preps

- + Minimally invasive Class II MO/DO cavity preps in just 30 to 60 seconds

Soft Tissue Indications

- + Incision, excision, vaporization, ablation and coagulation of oral soft tissues

*T&C Apply

Technical Specifications

Wavelength : Er;Cr:YSGG, 2780 nm

Power : 8 watts

Pulse Repetition Rates : 10-50 Hz

Pulse Energy : 300 mj



UNICORN DenMart 
Behind The Best Indian Dentists