

- ✦ Aluminium alloy chair to prevent rusting & reliability
- ✦ Synchronized movement of seat and backrest
- ✦ 6 Programs in the chair which include - 3 programmable working positions, preset auto return to zero, spitting and emergency switch
- ✦ Trendelenburg position for dental emergencies
- ✦ Independent up/down chair movements along with back/forth backrest movements
- ✦ Instrument control panel incorporates with 3 figure numerical display for greater data reading precision

LED LIGHT

- ✦ 3 axis rotation
- ✦ Inhouse temp. : 4,500k
- ✦ Light intensity : 30,000 Lux

CONTINENTAL DELIVERY UNIT

- ✦ 5 modules with balance arm hangers
- ✦ Separate feather-touch keypad
- ✦ Each instrument has a retractions system


V286

CE EN ISO 13485:2012


 Vector V286
Chair comes
with VMax 450
Suction System

- ✦ X-Ray viewer with automatic turn off
- ✦ Removable & autoclavable resting silicon pads
- ✦ Movable instrument tray below the delivery unit

WATER UNIT & ASSISTANT ARM

- ✦ High quality ceramic spittoon with programmable water input for spittoon and cup
- ✦ Assistant Instrument table, fitted on articulated arm for 2 suction and optional instruments
- ✦ 1 - 3 way syringe
- ✦ Dental stool with chair
- ✦ Spittoon bowl can be dismantled completely with removable cup & bowl fill unit

MULTI FUNCTIONAL FOOT CONTROL

- ✦ Joy stick for patient chair control
- ✦ Rinse & Zero position
- ✦ Light & cup fill control
- ✦ Chip air & chip water control
- ✦ Instrument speed & power control

