


- ✦ Aluminium casting chair
- ✦ New concept, aesthetic & modernized design
- ✦ Seamless & soft upholstery (Memory Foam) for patient comfort
- ✦ Synchronised movement of seat & backrest
- ✦ Chair can lift upto 250 kg weight
- ✦ Hydraulic delivery system from Japan
- ✦ Hand rest - with two direction rotation
- ✦ Safety with protection against any obstruction when chair goes down
- ✦ Double articulated headrest with pull & push button for pediatric & handicapped patients
- ✦ All connection in the base of chair
- ✦ 2 customizable & 1 last working position Program
- ✦ Anti retraction valve



The seat can be rotated about 30 degree, More convenient for maintenance & multiple operations.



DELIVERY UNIT


- ✦ Delivery unit with 3 air rotor point & one point for 3-way syringe
- ✦ All the control of chair are on unit with touch panel

WATER UNIT

- ✦ 90° Movements best for peado & old age patients
- ✦ Removable cuspidor, easy to be disinfected
- ✦ Air & water suction filter to prevent drainage block

ASSISTANT SIDE

- ✦ Three points - one for 3 way syringe & 2 for low high vacuum suction
- ✦ Chair control on assistant side.




DENTIST STOOL

- ✦ With adjustable up & down position & also backrest can be adjusted upto 250 degree

DCI LIGHT

- ✦ 3 axis rotation (For easy positioning and projection)
- ✦ 3-intensity : 16000 Lux, 26000 Lux & 38000 Lux
- ✦ Precise illumination Pattern : 3" X 6" (Prevents Patient Glare)
- ✦ Color Temperature : 5000K (Accurate shade matching)
- ✦ Color Rendering Index (CRI) : 90
- ✦ Shadow free light (Less eye fatigue & Increase Visibility)
- ✦ Reflective technology : Uses 500 reflective facet (Uniform light pattern)
- ✦ Durable Design : Heavy-duty structure & handle design
- ✦ Hermetically sealed front screen makes external cleaning easier and stops any dust or spring fluid getting inside the light



The seat can be rotated about 30 degree, More convenient for maintenance & multiple operations.



NEW

