



## ENDODONTIC EQUIPMENT

### ENDO-MATE AT

STYLISH & COMPACT ENDO MICROMOTOR

- ✦ Fits all major brands of Ni-Ti Rotary Files
- ✦ Torque control, with integrated auto-reverse function
- ✦ Compact Stylish Control/ Monitor pad
- ✦ Large screen Size for immediate visibility
- ✦ Flag control panel for simple & user friendly operation
- ✦ Light weight, Smart hand piece offering excellent visibility
- ✦ Convenient On/Off Switch
- ✦ Power supply AC 230 V 50/60 Hz
- ✦ **Nos of programs 9**
- ✦ **Max. Torque 7 Ncm (with 20:1 head)**
- ✦ **Speed range 100-13,000 min<sup>-1</sup>**
- ✦ **Dimensions W268 x D230 x H130 (mm)**
- ✦ **Weight 342 gms only (Including motor handpiece & cord)**
- ✦ **AC operated**

Endodontics has never  
been so **EASY...!!**

