



ENDODONTIC EQUIPMENT

ENDO-MATE DT

ENDO MICROMOTOR WITH TORQUE CONTROL AND AUTO REVERSE

- ✦ Variable Speed
- ✦ 9 programmable memory
- ✦ 6 head angles for ease of treatment
- ✦ Push type Ultra miniature head
- ✦ On/Off switch on hand piece
- ✦ Wide screen provides high visibility
- ✦ Compact desktop unit for all types of rotary files
- ✦ Speed range : 100-13,000 min⁻¹
(With 20:1, 4:1, 1:1, 16:1 geared heads)
- ✦ Dimensions : W92 x D148 x H124 (mm)
- ✦ Weight : Unit 456g/Motor handpiece & cord 92g
- ✦ Torque Level : 0.8~6.5 Ncm (When using 20:1 Handpiece)
- ✦ Battery / AC operated



EC ISO 9001:2008 ISO 13485:2003 ISO 14001:2004