

- ✦ Double articulated headrest for handicapped & pediatric patients, headlock system with lever

SHADOWLESS LED DENTAL LIGHT


- ✦ Intensity : 8000 - 35000 Lux
- ✦ Head made of resistant material, with a rotation of 620 degree, propitiating lightness, durability and ample mobility in all positions.
- ✦ Multifaceted mirror with multicoating treatment, generating innumerable light sources

MULTIFUNCTION FOOT CONTROL

- ✦ Four work programs and memory with the status of the dental light
- ✦ Spitting position and return to the last position
- ✦ Activation and adjustment of intensity of the dental light
- ✦ Automatic return to zero position
- ✦ Rise and lowering of the seat and backrest
- ✦ Movements blocking of dental chair with emergency light indication

OVERHEAD DELIVERY UNIT

- ✦ With pneumatic locking arm to prevent any unnecessary displacement
- ✦ Feather touch control panel for easy touch play operations
- ✦ Allow up to five working terminal
- ✦ Three point for air rotor
- ✦ One 3 way syringe
- ✦ Removable & autoclavable silicon pads for complete germ free operations
- ✦ Pressure gauze on the side of the delivery unit
- ✦ With inbuilt LED X-Ray Viewer



WATER UNIT

- ✦ 60 degree movable water unit
- ✦ Removable spittoon for easy cleaning
- ✦ Cup filler with electrical activation
- ✦ Proximity Sensor: Auto control the spittoon water
- ✦ Two points one for Low vacuum suction and other for high vacuum suction cannulas

CE ISO 13485:2012 + AC:2012 / ISO 9001:2008

ZERO BACK ACHE STOOL



- ✦ Available with all Gnatus Chairs. Helps prevent Spondylitis & low back pain. GNATUS has developed with considerable in R&D investment
- ✦ Up & down soft movements
- ✦ Adjustments of height & tilt are made through a side lever
- ✦ Adjustable from 50 to 70 cm of height with an adjustable ring
- ✦ The backrest with tilt & height adjustment which provides comfort to the lumbar area